

LED BULB FOR CLAR MIROR AND CHARDON MIROR

Ref. 570155



Powerful and long-lasting lighting for your examinations with Spengler Clar and Chardon mirrors.

Description

LED bulb for clar mirror and SPENGLER thistle mirror

31/10/2025 21:48:01 – For information purposes only, can be subject to change

LED BULB FOR CLAR MIROR AND CHARDON MIROR

Ref. 570155



Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446000073	HS Code	90189084
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	1	Packed unit / Height (cm)	1
Packed unit / Length (cm)	2	Packed unit / Weight	0.002

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---