

REPLACEMENT HEAD BLUEBERRY FOR OTOSCOPE SMARTLED 5500



Ref. 570518B

A reliable replacement head to extend the life of your Smartled 5500 otoscope.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446041779	HS Code	90189042
Country	Pakistan		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.046		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		