

RECHARGING DOCK FOR SMARTLED 5500-R



Ref. 579534

Fast and convenient charging to keep your Smartled otoscope always ready for use.

Description

This charging base is designed for the Spengler Smartled 5500-R otoscope. It provides a full charge in 2.5 hours via the included power adapter and can also be connected to a computer using a standard USB-C cable (not included). Comes with a power adapter and interchangeable plugs (European, American, British, and Chinese). 3-year warranty.

Features

| | | | |
|------------|---------------|----------------|----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3700446032913 | HS Code | 85366990 |
| Country | Pakistan | | |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|-------|---------------------------|-------|
| Packed unit / Width (cm) | 14 | Packed unit / Height (cm) | 8 |
| Packed unit / Length (cm) | 17 | Packed unit / Weight | 0.354 |
| Bare Unit / Weight (kg) | 0.354 | | |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|