

**CABLE USB-C - USB MALE 2M POUR  
SMARTLED 5500-R****Ref. 579535**

A sturdy and convenient cable for easy connection to the SMARTLED 5500-R.

## Description

This 2-meter USB-C to USB male cable ensures fast and reliable connectivity with the SMARTLED 5500-R monitor. Its extended length provides flexibility in setup while maintaining stable data transmission.

## Features

|            |               |                |            |
|------------|---------------|----------------|------------|
| CE Marking | NONE          | Medical device | No         |
| Bar code   | 3661040216782 | HS Code        | 9018110000 |
| Country    | China         |                |            |

### DIMENSIONS AND WEIGHT

|                           |     |                           |     |
|---------------------------|-----|---------------------------|-----|
| Packed unit / Width (cm)  | 14  | Packed unit / Height (cm) | 8   |
| Packed unit / Length (cm) | 17  | Packed unit / Weight      | 354 |
| Bare Unit / Weight (kg)   | 354 |                           |     |

### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -20 | Max temperature  | 80 |
| Minimum humidity    | 5   | Maximum humidity | 90 |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|