

FLOATING DIAPHRAGM BIG HEAD CHEST PIECE FOR RAPPAPORT TYPE STETHO



Ref. 605930

Clear and precise auscultation with a diaphragm tailored for Rappaport chestpieces.

Description

Spare part for Rappaport type stethoscopes.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446034757	HS Code	90181990
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.001

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---