

FITTING RINGS FOR SMALL HEAD CHEST PIECE FOR RAPPAPORT TYPE STETHO



Ref. 605941

Precise and durable fastening for the small chestpiece of your Rappaport stethoscope type.

Description

Spare part for Rappaport type stethoscopes.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446034788	HS Code	90181910
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.006

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---