


## PROBE FOR VASCULAR DOPPLER / 5 MHZ

Ref. 743311



CE  
0123



A precise and ergonomic 5 MHz probe for vascular assessments with the Spengler Doppler.

### Description

5 MHz probe - vascular probe for doppler

31/10/2025 21:49:48 – For information purposes only, can be subject to change

## PROBE FOR VASCULAR DOPPLER / 5 MHZ

Ref. 743311



### Features

|                   |               |                |           |
|-------------------|---------------|----------------|-----------|
| CE Marking        | CE0123        | Medical device | Class IIa |
| Bar code          | 3700446008178 | HS Code        | 90181990  |
| Warranty (months) | 12            | Country        | China     |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.145 |

#### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature  | 55 |
| Minimum humidity    | 5   | Maximum humidity | 93 |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|