

POWER ADAPTER

Ref. [DDAPDB07](#)

Reliable power supply to keep your BM-100A bilirubinometer ready at all times.

Description

The power adapter for the BM-100A bilirubinometer (ref. DDAPDB07) enables charging through its dedicated base. It provides stable and safe power, essential for repeated measurements in neonatal care. Not classified as a medical device, CE marked.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040220949	HS Code	85444
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---