

BASE

Ref. [DDAPDB09](#)

A stable and discreet base for safely positioning the bili BM-100A.

Description

This base is designed to hold the bili BM-100A securely, ensuring a fixed and safe setup. It supports optimal positioning of the device during phototherapy sessions.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040220963	HS Code	9018909911
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---