

SOFT EARPIECE FOR IDEAL+ (PAIR), RED

Ref. DS8009_V



Soft large-size eartips for comfortable and hygienic auscultation with the IDEAL+.

Description

=====

31/10/2025 21:51:53 – For information purposes only, can be subject to change

SOFT EARPIECE FOR IDEAL+ (PAIR), RED

Ref. DS8009_V



Features

CE Marking	NONE	Medical device	No
Bar code	3661040073996	HS Code	90189012
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	1
Packed unit / Weight	0.004	Bare Unit / Weight (kg)	0.006

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Green	Packed unit / Quantity	1
--------	-------	------------------------	---