

**CT50 RECHARGEABLE LITHIUM ION
BATTERY****SUNTECH**Ref. DSU222210

Extended autonomy and reliable performance for your SUNTECH CT50 monitor.

Description

This lithium battery is specifically designed for the SUNTECH CT50 monitor. It provides long-lasting use, fast recharging, an optimal reliability for demanding medical environments.

Features

CE Marking	NONE	Medical device	No
Bar code	3661040214283	HS Code	850650
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	6	Packed unit / Height (cm)	3
Packed unit / Length (cm)	20	Packed unit / Weight	0.534

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---