

CABLE TANGO ECG INTERFACE FOR GE CASE, CASE 800

SUNTECH

Ref. DSU230742

A reliable interface cable to connect TANGO with your GE Case ECG systems.

Description

This ECG cable enables seamless integration between the TANGO monitor and GE Case or Case 800 systems. It ensures accurate cardiac data transmission for stress testing or ECG monitoring workflows.

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3700446012397	HS Code	9018110000
Warranty (months)	3	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	2
Packed unit / Length (cm)	20	Packed unit / Weight	0.15

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---