

ADULT SUNTECH ARMBAND

Ref. DSU230816

SUNTECH

A reliable adult cuff for comfortable and accurate blood pressure monitoring.

Description

This SunTech adult cuff is designed to provide accurate blood pressure readings with optimal fit and comfort. Compatible with various devices, it supports efficient use in clinical environments.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446012519	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.074

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---