

INFANT 8-13CM X20

Ref. DSU231080

SUNTECH



Double-tube infant cuffs for accurate and hygienic blood pressure monitoring in pediatric care.

Description

Set of 20 Suntech "SOFT" single-patient armbands, orange.
Double tubing, Small child size 8-13cm, screw connection.

31/10/2025 21:38:10 – For information purposes only, can be subject to change

**INFANT 8-13CM X20**Ref. DSU231080**SUNTECH**

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014490	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.653

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---