

CHILD 12-19CM X20

Ref. DSU231081

SUNTECH



CE



Single-patient pediatric cuffs for reliable and hygienic monitoring in clinical settings.

Description

Set of 20 Suntech "SOFT" single-patient armbands, green.
Double tubing, Child size 12-19cm, screw connection.

31/10/2025 21:42:08 – For information purposes only, can be subject to change

CHILD 12-19CM X20

Ref. DSU231081

SUNTECH

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014506	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.792
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		