

SUNTECH CUFF FOR SINGLE PATIENT NEONAT N°2 (4-8CM)

SUNTECH

Ref. [DSU231101](#)

Single-use neonatal cuffs for gentle and accurate monitoring from day one.

Description

This pack of 20 SunTech single-use cuffs, size No. 2 (4–8 cm), is designed for newborns. It ensures precise blood pressure measurement while minimizing cross-contamination risks, ideal for neonatal intensive care settings.

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935103421	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	16
Packed unit / Length (cm)	11	Packed unit / Weight	0.25

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---