

SET OF 20 SUNTECH MONOBLOC CUFFS SINGLE PATIENT NEONATAL N°5 (8-15CM)

SUNTECH

Ref. [DSU231104](#)

Single-use neonatal cuffs for safe and tailored blood pressure monitoring from day one.

Description

Designed for newborns, these single-use cuffs provide accurate blood pressure readings while minimizing cross-contamination risks. The one-piece design ensures quick and secure application.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014575	HS Code	9018905040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Packed unit / Height (cm)	11
Packed unit / Length (cm)	31	Packed unit / Weight	0.434

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---