

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE CHILD (12-19CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. [DSU231111](#)



CE



Durable pediatric cuffs for reliable monitoring with metal Bayonet connector.

Description

Reusable nylon cuff sold in boxes of 5 units.
Single tubing, bayonette connection 12-19cm.

31/10/2025 21:41:24 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC ARMBANDS REUSABLE CHILD (12-19CM) - BAYONET METAL CONNECTOR

SUNTECH

Ref. [DSU231111](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041008799	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	9	Packed unit / Height (cm)	14
Packed unit / Length (cm)	11	Packed unit / Weight	0.25

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---