

FLEXIBLE CHILD DIGITAL SPO2 SENSOR (90CM CABLE)

SUNTECH

Ref. [DSU250904](#)

Reliable and gentle oxygen monitoring for children, compatible with SunTech devices.

Description

Flexible child digital SpO2 sensor (90 cm cable)

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	8698870322296	HS Code	90181990
Country	Turkey		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	2	Packed unit / Height (cm)	16
Packed unit / Length (cm)	14	Packed unit / Weight	0.01

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---