

**AXILLARY ORAL TEMPERATURE PROBE
FOR COVIDIEN FILAC3000 TEMP MODULE****SUNTECH****Ref. DSU551463**

Fast and reliable oral or axillary temperature measurement with the Filac 3000 module.

Description

This temperature probe is compatible with the Covidien Filac 3000 module. Designed for oral or axillary measurements, it delivers fast and accurate readings, ideal for demanding clinical settings.

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	0120884521008134	HS Code	9018.19.5500
Warranty (months)	3	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	8	Packed unit / Height (cm)	5
Packed unit / Length (cm)	15	Packed unit / Weight	0.086

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---