

F3000 DISPOSABLE TEMP PROBE COVERS (25 BOXES/TRAY, 20 COVERS/BOX)

SUNTECH

Ref. [DSU551467](#)

Disposable hygienic protection for safe and compliant temperature measurements.

Description

25 boxes of 20

Features

CE Marking	CE0413	Medical device	Class I
Bar code	1088452116202	HS Code	9018.19.5500
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Packed unit / Height (cm)	10
Packed unit / Length (cm)	29	Packed unit / Weight	0.622

SPECIFIC STORAGE CONDITIONS

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---