

OPD-CHILD CUFF (5 PIECES 12-19CM) BRAVO MINI SUNTECH

SUNTECH

Ref. [DSU552555](#)

Pediatric cuffs designed for accurate blood pressure monitoring with the Bravo Mini ABPM.

Description

This set of 5 pediatric OPD cuffs is compatible with the Suntech Bravo Mini ABPM. Designed for arm sizes from 12 to 19 cm, they provide optimal comfort and reliable blood pressure readings for young patients during ambulatory monitoring.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040217707	HS Code	90189010
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---