

POWER CORD EUROPE PRO VS CHECK

SUNTECH

Ref. DSU552760

Safely powers your VS Check devices with European mains compatibility.

Description

Spare part for tourniquets

Features

CE Marking	CE0413	Medical device	Class I
Bar code	3700446033590	HS Code	85049099
Warranty (months)	3	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	8	Packed unit / Height (cm)	4
Packed unit / Length (cm)	222	Packed unit / Weight	0.27

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---