

Ref. MONI011D



CE



30/10/2025 20:41:27 – For information purposes only, can be subject to change

=====

Ref. MONI011D

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002483	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.06
Bare Unit / Weight (kg)	0.055		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES

📅 30/10/2025 20:41:27 – For information purposes only, can be subject to change