

Ref. MONI011G



CE



30/10/2025 21:00:15 – For information purposes only, can be subject to change

=====

Ref. MONI011G

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
HS Code	9018193090	Warranty (months)	6
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085
Bare Unit / Weight (kg)	0.08		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES

 30/10/2025 21:00:15 – For information purposes only, can be subject to change