

Ref. MONI011



CE



31/10/2025 21:37:56 – For information purposes only, can be subject to change

Ref. MONI011

Description

Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002421	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.08
Bare Unit / Weight (kg)	0.075		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES

31/10/2025 21:37:56 – For information purposes only, can be subject to change