

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513210



CE
0459



Light and robust, the blood pressure monitor that combines style and performance.

Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

03/11/2025 23:26:08 – For information purposes only, can be subject to change

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513210



Features

CE Marking	CE0459	Medical device	No
Bar code	3700446009953	HS Code	90189010
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	7
Packed unit / Length (cm)	19	Packed unit / Weight	0.5
Bare Unit / Weight (kg)	0.4		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--