

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513239



CE
0459



Light and robust, the blood pressure monitor that combines style and performance.

Description

MANOMETRE EN CONTROLE QUALITE : VERIFIEZ QUE VOTRE LOT EST ACCEPTE, AVERTIR POUR CONTROLE PAR ECHANTILLONNAGE

03/11/2025 22:50:53 – For information purposes only, can be subject to change

LIAN NANO/SILVER RING/ BLACK NYLON / MULTI CUFFS (L, M, S)

Ref. 513239



Features

CE Marking	CE0459	Medical device	No
Bar code	3700446024703	HS Code	90189010
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.498

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---