

HYDROCELE TROCART 3.5 MM, 19CM

Ref. IT11411



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031682	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	6	Packed unit / Height (cm)	3
Packed unit / Length (cm)	22	Packed unit / Weight	0.039
Bare Unit / Weight (kg)	0.039		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--