

=====

Ref. GCDMS105

=====



Description

=====

Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | CE            | Medical device            | Class I    |
| Bar code                    | 3700468208143 | HS Code                   | 9018908400 |
| Warranty (months)           | 24            | Country                   | France     |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | nc         |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |