

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 163 AND CLAMP/ADULT (M)

Ref. 337883



Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Datascope monitor

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446000172	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.125
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		