



Ref. TS53000S



CE
0459



07/11/2025 19:39:26 – For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

Ref. TS53000S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215486	HS Code	90189010
Country	France		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	25
Packed unit / Length (cm)	16	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

07/11/2025 19:39:26 – For information purposes only, can be subject to change