

Ref. MONI010D

 Spengler



CE
0123



Description

 07/11/2025 19:44:07 – For information purposes only, can be subject to change

Ref. **MONI010D**



Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002537	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
Bare Unit / Weight (kg)	0.05		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES

07/11/2025 19:44:07 – For information purposes only, can be subject to change