

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT

Ref. BCASPI38

CA-MI



CE
0123

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

Intended for the use of professionals.

Description

C30 Suction Unit, Flovac® single-use pouches, 2 x 2000 ml, Foot

10/11/2025 21:06:10 – For information purposes only, can be subject to change

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT

CA-MIRef. BCASPI38

control-Foot

control

Surgical Suction Unit on wheels and mobile, intended to drain body fluids and to practice oral, nasal and tracheal aspiration in adults and children.

Intended for professional use and equipped with a powerful suction system that can drain up to 40 l/min max, Askir C30 is equipped with two collection jars and a system on five casters making it easy to transport.

- Main system equipped with a pressure gauge (bar & kPa) and a vacuum regulator
- Available in two versions: Basic or FS
- Accessories for the Basic version:

-

Fluid collection jar with an anti-reflux valve

-

8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm

-

Antibacterial & hydrophobic filter (single patient)-

Power cord with Schuko plug

- Accessories for the FS version

:-

Fluid collection jar with an anti-reflux valve

-

8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm

-

Antibacterial & hydrophobic filter (single patient)

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT



Ref. BCASPI38

- Power

cord with Schuko plug

- Pedal

- Technical characteristics

:-

Engine: Oil-free piston pump requiring no maintenance

-

Power supply: 230V-50Hz

- Maximum vacuum (adjustable): -0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 40 l/min

-

Noise level: 61 dBA

-

Weight: 6.2 kg-

Dimensions: 32 x 99 x 30 cm

SUCTION UNIT CA-MI ASKIR C30 FS BOTTLE/ BAG UU FLOVAC 2X2L CDE/FOOT



Ref. BCASPI38

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139241	HS Code	90189084
Warranty (months)	24	Country	Italy

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	9.15
Bare Unit / Weight (kg)	6.5		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85

OTHER FEATURES