

DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN

Ref. SP1644P



CE
0123



DISPLAY OF 12 THERMOMETER TEMPO 10 FLEX GREEN



Ref. SP1644P

Display of 12 Thermometer Tempo 10 Flex green

Description

Display of 12 thermometers TEMPO 10 FLEX

The thermometer TEMPO 10 FLEX is a thermometer indicating the body temperature in 10 seconds and having a flexible tip

- Fever alarm
- Automatic shutdown
- Water resistant
- Last measure memory
- 2 colors available purple or green
- Comes with: storage case, button battery (1.5V).

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446027537	HS Code	9025190090
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	11	Packed unit / Height (cm)	15
Packed unit / Length (cm)	14	Packed unit / Weight	0.624
Bare Unit / Weight (kg)	0.12		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--