

Ref. **BM104S**

Description

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468215752	HS Code	9018901000
Warranty (months)	12		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES