

DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. GBDM102

DDM



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE



13/11/2025 21:52:02 – For information purposes only, can be subject to change

DECONTAMINABLE ARM CUFF MONOBLOC STRAIGHT NEW BORN TOURNIQUET DDM

Ref. GBDM102



Monobloc straight arm cuff with one bladder and a ribbon

Description

Monobloc PU décontaminable arm cuff are available in straight design. Made of a single monobloc unit with one bladder and ribbon to secure the Velcro attachment.

They are available in a range of different sizes easily indentifiable by color coding of the tubing. (A : bladder width, B: bladder length, C: total length)

		A	B	C
• New born arm cuff	White tube	3.5	29	38
• Child arm cuff	Pink tube	4.5	35	48
• Small adult arm cuff	Parme tube	6	46	61
• Adult arm cuff	Sky-blue tube	8	46	67
• Large adult arm cuff	Turquoise tube	10	47	70

Monobloc arm cuffs can not be sterilized. They can be dry cleaned or immersed paying attention to maintain the tube connector out of water.

**DECONTAMINABLE ARM CUFF MONOBLOC
STRAIGHT NEW BORN TOURNIQUET DDM**

Ref. GBDM102

**Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468207917	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES