

TROLLEY WITH WHEELS ELITE BAGS

Ref. MECHAR02

 **ELITE BAGS**



 13/11/2025 21:17:43 – For information purposes only, can be subject to change

TROLLEY WITH WHEELS ELITE BAGS

Ref. MECHAR02



External trolley, foldable, robust and lightweight.

Description

It folds quickly thanks to a button located on the handle. It has an elastic rope that allows you to completely fix the load.

Materials: steel tubes, ABS handle, PP base and PVC wheels

Wheel dimensions: Ø 77 x 22 mm

Color: black

Weight: 1.35 kg- Outdoor measures: 23 x

96.5 x 32 cm (folded: 23 x 44 cm)

Recommended max. capacity: 30 kg

- Can be used with models: Call, Podia, Doctor, Sport, Medic, Cool, Practi, Vet, Urban and Doctor

TROLLEY WITH WHEELS ELITE BAGS

Ref. MECHAR02



Features

CE Marking	NONE	Medical device	No
Bar code	8436564540621	HS Code	4202129900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	25	Bare Unit / Width (cm)	23
Packed unit / Height (cm)	11	Bare Unit / Height (cm)	32
Packed unit / Length (cm)	48	Bare Unit / Length (cm)	96.5
Packed unit / Weight	1.8	Bare Unit / Weight (kg)	1.35

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---