



Ref. TS50004S



CE
0459



14/11/2025 20:27:37 – For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

Ref. TS50004S



Description

Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215455	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	20
Packed unit / Length (cm)	13	Packed unit / Weight	0.485
Bare Unit / Weight (kg)	0.325		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

OTHER FEATURES

14/11/2025 20:27:37 – For information purposes only, can be subject to change