

LIAN METAL HAND-HELD MANOMETER

Ref. 310045



Description

Dials 0 to 300 mmHg

17/11/2025 21:08:58 – For information purposes only, can be subject to change

LIAN METAL HAND-HELD MANOMETER

Ref. 310045



Features

CE Marking	NONE	Medical device	Class Im (measuring function)
Bar code	3700446000141	HS Code	90189010
Warranty (months)	60	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.38
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		