

CUFF FOR NEOTENS® ELECTRONIC BLOOD PRESSURE MONITOR - SPENGLER

Ref. 527052



CE



17/11/2025 21:06:24 – For information purposes only, can be subject to change

CUFF FOR NEOTENS® ELECTRONIC BLOOD PRESSURE MONITOR - SPENGLER

Ref. 527052

Range of cuffs for NEOTENS® electronic blood pressure monitor

Description

NEOTENS® blood pressure cuffs are very comfortable and come in three sizes:

- Adult cuff M/L (22-42 cm) for NEOTENS® - Reference: 527051
- Child cuff S (16-24 cm) for NEOTENS® - Reference: 527050
- Adult cuff XL (40-55 cm) for NEOTENS® - Reference: 527052

Features

CE Marking	CE	Medical device	Class I
Bar code	6976054540158	HS Code	90262040
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-25	Max temperature	60
Minimum humidity	NA	Maximum humidity	90

OTHER FEATURES