


VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)

Ref. 512200



CE
0459



23/11/2025 11:07:40 – For information purposes only, can be subject to change

**VAQUEZ LAUBRY NANO / VELCRO CUFF /
BLACK / CHILD (S)****Ref. 512200**

=====

Description

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007263	HS Code	90189010
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	20	Packed unit / Weight	0.8

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---