

## BABY TO TODDLER ELECTRONIC SCALE

Ref. EHANDI3



23/11/2025 10:52:58 – For information purposes only, can be subject to change

## BABY TO TODDLER ELECTRONIC SCALE

Ref. EHANDI3



**4-in-1 baby and child scales: ideal for paediatricians and general practitioners**

### Description

Baby and child scale that allows to weigh and measure babies with a removable tray or to weigh and measure young children from 30 to 160 cm without the tray. •  
 3 weighing modes: classic mode, milk absorption and mom + baby weight • Large capacity up to 100 kg: to weigh mom + baby •  
 Large ergonomic weighing tray area resistant and optimized to weigh infants (600(L) x 326(W) x 91(H) mm) •  
 4 touch keys including manual hold function  
 • Hygienic, Easy to clean, stable with 4 non-slip rubber feet, easy to carry (Carrying bag sold separately) •  
 Conversion of units possible to kg/lb and cm/ft •  
 Large LCD display (79 x 28 mm) with blue backlight, easy to read and automatically rotates 180° when weight is greater than 20 kg for convenient reading  
 • Power: 3 x 1.5V AAA batteries (not included)  
 • Auto On-Off function for power saving

### Features

CE Marking	CE	Medical device	No
Bar code	3661040216669	HS Code	84231090
Warranty (months)	12	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	34	Packed unit / Height (cm)	16
Packed unit / Length (cm)	63	Packed unit / Weight	1.8
Bare Unit / Weight (kg)	1.6		

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--