

LITESTICK PENLIGHT SPENGLER

Ref. 570621



LITESTICK PENLIGHT SPENGLER

Ref. 570621



Pen-shaped medical light, designed for close-up clinical examinations.

Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

Features

CE Marking	CE	Medical device	No
Bar code	3700446009861	HS Code	9018908400
Warranty (months)	24		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.015		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	40
---------------------	----	-----------------	----

OTHER FEATURES

Colors	Blue	Packed unit / Quantity	1
--------	------	------------------------	---