

# CABLE PATIENT ECG TANGO

Ref. [DSU230716](#)

SUNTECH

Cable patient ECG Tango

## Description

ECG cable for Tango® automatic monitor

## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE0413        | Medical device | Class IIa  |
| Bar code          | 3700446012298 | HS Code        | 9018110000 |
| Warranty (months) | 3             | Country        | China      |

### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc    |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | 0.217 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|