

## EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

Ref. [MEEMER02\\_B](#)

**EMERGENCY'S™**



24/11/2025 21:43:38 – For information purposes only, can be subject to change

## EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

**EMERGENCY'S™**

Ref. MEEMER02\_B

**Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment**

### Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures

-

Reflective strips- Washable case

- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg

-

Available in waterproof version

- Sold empty

- Available in blue (MEEMER02\_B), yellow (MEEMER02\_J), black (MEEMER02\_N), orange (MEEMER02\_O), pink (MEEMER02\_RC) and red (MEEMER02\_R).

## EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

**EMERGENCY'S™**

Ref. MEEMER02\_B

### Features

|                             |               |                         |           |
|-----------------------------|---------------|-------------------------|-----------|
| CE Marking                  | NONE          | Medical device          | No        |
| Bar code                    | 8436564540874 | HS Code                 | 420299000 |
| Warranty (months)           | 24            | Country                 | China     |
| DIMENSIONS AND WEIGHT       |               |                         |           |
| Packed unit / Width (cm)    | 24            | Bare Unit / Width (cm)  | 25        |
| Packed unit / Height (cm)   | 22            | Bare Unit / Height (cm) | 27        |
| Packed unit / Length (cm)   | 30            | Bare Unit / Length (cm) | 44        |
| Packed unit / Weight        | 1.3           | Bare Unit / Weight (kg) | 1.29      |
| SPECIFIC STORAGE CONDITIONS |               |                         |           |
| Minimum temperature         | NA            | Max temperature         | NA        |
| Minimum humidity            | NA            | Maximum humidity        | NA        |
| OTHER FEATURES              |               |                         |           |
| Colors                      | Blue          | Packed unit / Quantity  | 1         |