

## MALE TUNIC, BOUGAINVILLEA LYOCELL , S0, HOLTEX

Ref. WMBOU07\_90



📅 24/11/2025 21:20:20 – For information purposes only, can be subject to change

## MALE TUNIC, BOUGAINVILLEA LYOCELL , S0, HOLTEX

Ref. WMBOU07\_90



### The eco-responsible men's medical tunic

#### Description

- White men's tunic with short sleeves,
- Eco-responsible: 50% Lyocell (eucalyptus fiber) and 50% Polyester,
- Officer collar,
- Twill, 215 g / m2,
- 3 plated pockets,
- Buttoning by pressure.
- Available in :
  - Size 0 (WMBOU07\_90)
  - Size 1 (WMBOU07\_91)
  - Size 2 (WMBOU07\_92)
  - Size 3 (WMBOU07\_93)
  - Size 4 (WMBOU07\_94)
  - Size 5 (WMBOU07\_95)

The Holtex eco-friendly and responsible medical tunic brings extra comfort and strength to a classic polycotton tunic: superior quality, increased strength, easy maintenance and maximum ease. Its cut is traditional and modern at the same time. Its officer collar, its short sleeves mounted lapels and its 3 plated pockets are nicely finished with a white border which makes this men's tunic a very neat garment.

Lyocell is a natural material, made from eucalyptus wood pulp.

- Soft, lightweight, comfortable and slightly elastic,
  - Easy to maintain, wrinkle resistant, it does not lint.
  - Antibacterial, hypoallergenic, breathable, it absorbs moisture (up to 50% more than cotton): it is ideal for sensitive skin.
  - Environmentally friendly: Eucalyptus cultivation requires much less water and pesticides than cotton.

**MALE TUNIC, BOUGAINVILLEA LYOCELL ,  
S0, HOLTEX**

Ref. WMBOU07\_90



## Features

|                             |                      |                           |              |
|-----------------------------|----------------------|---------------------------|--------------|
| CE Marking                  | NONE                 | Medical device            | No           |
| Bar code                    | 3661041009475        | HS Code                   | 6211433100   |
| Country                     | Tunisia              |                           |              |
| DIMENSIONS AND WEIGHT       |                      |                           |              |
| Packed unit / Width (cm)    | nc                   | Packed unit / Height (cm) | nc           |
| Packed unit / Length (cm)   | nc                   | Packed unit / Weight      | nc           |
| SPECIFIC STORAGE CONDITIONS |                      |                           |              |
| Minimum temperature         | NA                   | Max temperature           | NA           |
| Minimum humidity            | NA                   | Maximum humidity          | NA           |
| OTHER FEATURES              |                      |                           |              |
| Bleaching                   | Bleaching prohibited | Washing                   | Wash at 60°C |
| Model                       | Male model           | Cleaning                  | Dry cleaning |
| Packed unit / Quantity      | 1                    | Size                      | Size 0       |