

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. WSAM\_47BF

medi Plogs



## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. [WSAM\\_47BF](#)



### Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

### Features

|            |               |                |          |
|------------|---------------|----------------|----------|
| CE Marking | NONE          | Medical device | No       |
| Bar code   | 3661040100425 | HS Code        | 64029993 |
| Country    | Germany       |                |          |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |      |
|---------------------------|----|---------------------------|------|
| Packed unit / Width (cm)  | 22 | Packed unit / Height (cm) | 12   |
| Packed unit / Length (cm) | 31 | Packed unit / Weight      | 1.15 |
| Bare Unit / Weight (kg)   | 1  |                           |      |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|        |           |           |    |
|--------|-----------|-----------|----|
| Colors | Navy blue | Shoe Size | 47 |
|--------|-----------|-----------|----|