

REFLEX HAMMER ACCORDING TO TRÖMNER, ALUMINIUM - LUXAMED

Ref. M1.111.112

LUXAMED®



CE



REFLEX HAMMER ACCORDING TO TRÖMNER, ALUMINIUM - LUXAMED

LUXAMED®

Ref. M1.111.112

Reflex Hammer according to Trömnér | large

Description

- Reflex hammer made of aluminium / brass
 - Optimal weight distribution
 - With neurological sensor tip
 - Length: 23 cm | weight: 245 g
 - Without brush due to hygienic reasons
- Available in white (M1.111.915) and black (M1.111.112)

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960287	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	Black
--------	-------