

Ref. E1.314.012

**LUXAMED®**



CE



## Description

28/11/2025 19:48:13 – For information purposes only, can be subject to change

Ref. E1.314.012

LUXAMED®

## Features

|                             |                |                           |          |
|-----------------------------|----------------|---------------------------|----------|
| CE Marking                  | CE             | Medical device            | Class I  |
| Bar code                    | 04260381960416 | HS Code                   | 90189084 |
| Warranty (months)           | 24             | Country                   | Germany  |
| DIMENSIONS AND WEIGHT       |                |                           |          |
| Packed unit / Width (cm)    | nc             | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)   | nc             | Packed unit / Weight      | nc       |
| SPECIFIC STORAGE CONDITIONS |                |                           |          |
| Minimum temperature         | 5              | Max temperature           | 40       |
| Minimum humidity            | 10             | Maximum humidity          | 90       |
| OTHER FEATURES              |                |                           |          |
| Size                        | Size 4         |                           |          |